

# P2 filter for Respira Compact half mask

COD ARTICOL

SM0704017799999



## DESCRIZIONE

Possibility of filter sterilization up to 5 times. Soak in 60% alcohol, or by 80°C temperature for 20 minutes, or by ozone, or by UV. Nanofilters capture up to 50 nanometers, while the smallest virus has 200 nanometers. Much longer lifespan, thanks to the nanomembrane and one inhaling and exhaling valve around 4 times longer than usual half masks and full face masks filters. Package contains 3 pcs.

## STANDARDI

EN 143:2000/A1:2  
006

P2R

# P2 filter for Respira Compact half mask

COD ARTICOL

SM0704017799999

STANDARDE

EN 143:2000/A1:2006 / P2R

## CARACTERISTICI

Fără latex  
Fabricat in EU  
Fără silicon  
Greutatea articolului: 40g

Copyright © CERVA GROUP a.s .. Toate drepturile rezervate. Orice erori de imprimare, modificări de imagine și specificații ale produsului sunt rezervate.

Data modificării: 27 feb. 2023,  
16:36:05

CERVA GROUP a.s.,  
Průmyslová 483,  
252 61 Jeneč,  
Republica Cehă